


Dressing band cut to fit the wound
ADF HYDROCOLLOID DRESSIN



wound protection/ healing/ prevention of infection



standard	10ea, 4in*4in(10*10cm),
use	Peel off the tape according to the wound and attach the adhesive surface evenly
effect	Absorption and protection of secretions such as wounds and affected areas
ingredient	Hydrocolloid sheet
Storage method	airtight container, Room temperature storage(1~30°C)



Hydrocolloid protects wounds from water, dust, and bacteria and absorbs exudates from wounds to create a wet environment. It helps recover wounds by maintaining proper humidity of wounds.



- ① **Waterproof function**
- ② **smooth air permeation of the wound**
- ③ **Bacterial Blocking**
- ④ **Creating a Wet Environment for Wound Healing**

Homepage : (KR) <https://www.reacnf.com> ,(ENG) <https://kr.kompass.com/c/reacnf-co-ltd/kr033161/>

Video : <https://youtu.be/1gLLM-Sy5qw/>