







DermaQueen

Anti-aging Injection

DermaQueen can inject drug into the skin easily and quickly with a multi-needle consisting of several needles while holding the skin with the negative pressure during injection.

Also, it can evenly inject drugs into the desired skin layer during the procedure, increasing effectiveness with sophisticated procedures and enabling continuous drug injection, which lead to the reduction of drug loss and pain during the procedure.

Features

-  All types of syringes, including 10cc sized syringes, can be installed
-  All function controls implemented digitally
-  Adjustable injection length and quantity
-  Minimize drug loss

Gauge : 31G , 32G
Pin : 5 pin, 9 pin, 16 pin
Coming soon



Multi Needle







Portable type



Stand type

Surgical Blade (CE)

Features

-  Ensure the best quality with special heat-treatment
-  Sterilized by gamma ray and EO gas
-  Grinded keenly with newest auto-grinding machine
-  Made of special stainless steel





Use Handle No.3


Use Handle No.4



Chisel blade



Surgical Suture




SILK

Color

Size USP 8/0 ~ 1 | EP 0.4 ~ 4

Suture made of high-purity natural silk consist of 99% protein through degumming method, which is finely braided thin threads from the center at a low speed for high tensile strength.



NYLON

Color

Size USP 10/0 ~ 2 | EP 0.2 ~ 5

Nylon is made up of monofilaments and can minimize bacterial growth. It is made up of a single strand and has a smooth surface, making it easy to pass through the tissue.




POLYESTER

Color

Size USP 6/0 ~ 2 | EP 0.7 ~ 5

Braid threads made of polyester and coat them with silicone. Multi-filament provides excellent tensile strength and handling and knot stability. Gently passes through the tissue with the coated surface.



POLYPROPYLENE

Color

Size USP 8/0 ~ 1 | EP 0.4 ~ 4

Provides strong and stable support for tissue while providing continuous tensile strength to infected areas. Also, it is heat-resistant and has a low risk of infection. The adhesion to skin tissue is weak, so it is easy to remove the threads.

Suture for buried suture technique

Atraumatic Special Needle for non-incisional Blepharoplasty

Monofilament Nylon, 7/0, 40cm, Round 3/8 double
18mm, 20mm, 24mm / Black, Blue, White

DQ Cosmetic

Features



Fundamentally prevents loss of moisture to the external environment



Designed to relieve pain and prevent inconvenience in daily life



Silicone technology for exfoliating and moisturizing

Ingredients








Dimethicone/Vinyl Dimethicone Cross polymer
Dimethicone
Stearyl Dimethicone
Beeswax



•80ml

Johns Screw Needle

Features

-  The world's first needle length adjustable type (Adjustable range 0 ~ 4mm)
-  Ideal for vertical injection method for target dermis layer shooting
-  Korea patented / PCT pending application
-  Mainly used for injection of HA filler and Meso filler
-  Suitable for all types of syringes (up to 10ml size)



32G, 34G > **Coming soon**



White
27G x 7mm
(adjustable range 0 ~ 4mm)
100 Needles



Red
29G x 7mm
(adjustable range 0 ~ 4mm)
100 Needles



Pink
31G x 7mm
(adjustable range 0 ~ 4mm)
100 Needles

Johns Depressor

Features

-  Johns Depressor evenly distributes drug after the injection and allows uniform absorption

	Johns Injection Method	Linear Injection Method
Bleeding	Less	More
Pain	Painless	Painly
Tissue Damage	Small	More
Ecchymosis	Less	More
Lump	Less	More
Needle/Technique (Gauge, Length)	- Vertical Injection method (Evenly) - Exact Skin Thickness - Easy to handle	- Horizontal Injection method - Irregular Skin Thickness - Difficult to handle - Long length needle



2.5 mm cap
Full face



4 mm cap
Nasal tip
Augmentation



9 mm cap
Lip Augmentation



Full face



Full face



Malar



Chin